

User Acceptance Testing Guide

Role: Administrator — demo persona “Admin AIIA”

Document date: 12 September 2026

Application under test: **<https://ayusphere-three.vercel.app>**

Environment: Production demo (seeded demo data — no real patient data)

YOUR TEST CREDENTIALS

URL

Email

Password

This guide walks you, step by step, through every feature available to the **Administrator** role. No technical knowledge is required — follow the numbered steps, compare what you see with the screenshot, tick the checklist, and record Pass / Fail in the sign-off table at the end. All screenshots were captured from the live application.

AYUSH — Traditional Wisdom, Modern Evidence

About this role

What the Administrator can do: everything. The admin role holds every capability in the RBAC matrix — all PI powers (trials, participants, CRF entry *and approval*), the Ethics IEC gate, the PV adverse-event review walk, monitoring, alerts, copilot, documents, exports — *plus* the two admin-only areas: **Settings** (AI model configuration and user management) and the full role picture of the **Audit Trail** (shared with the Regulator).

This guide focuses on the admin-specific surfaces (Settings, Audit) plus a sweep of the shared ones; for deep tests of trials, CRF, ethics or PV flows, reuse the corresponding role guides — every step in them also works when logged in as admin.

What you will test

1. Logging in
2. Navigation overview & dashboard
3. Clinical Trials — full oversight
4. Settings — AI model configuration & user management
5. Audit Trail
6. Messages (internal communication)
7. Signing out

Before you start

General notes for testers:

- Use a modern desktop browser (Chrome or Edge recommended) at a comfortable window size.
- This is a shared demo environment — other testers may add data while you test; counts in your screenshots may differ slightly from the figures in this guide. That is normal and not a defect.
- Every action you take is recorded in the immutable audit trail — this is by design (ALCOA+ compliance).
- AI features (Doctor Note extraction and the AI Copilot) are powered by a configurable AI provider — currently Google Gemini (gemini-2.5-flash); an administrator can switch to Anthropic Claude under Settings. AI output can vary between runs; judge behaviour, not exact wording.
- If a page fails to load, refresh once before recording a failure (cold serverless starts can take a few seconds).

1 Logging in

What this is for. AyuSphere is protected by role-based access control (RBAC). Signing in with your role account determines exactly which pages and actions you can see — nothing more, nothing less.

STEPS

1. Open a browser and go to `https://ayusphere-three.vercel.app`. You are redirected to the login page.
2. In the **Email** field type `admin@aiaa.demo`.
3. In the **Password** field type `Demo@1234`.
4. Click the **Sign in** button.

Expected result. You are taken to the **Dashboard**. The top-right corner shows “Admin AIIA — Administrator”.

ACCEPTANCE CHECKLIST

- ☐ The login page loads with AyuSphere branding.
- ☐ A wrong password shows an error and does NOT log you in (try once, then use the correct one).
- ☐ The correct credentials open the dashboard.
- ☐ Your name and role are shown in the top-right corner.




Fig. 1. The AyuSphere login page. Captured from the live system on 12 September 2026.

2 Navigation overview & dashboard

What this is for. As admin the sidebar shows the entire platform — the union of every role's menu — which is itself a key RBAC test.

STEPS

1. After logging in, check the left sidebar. As Administrator you should see all 17 items: **Dashboard, Messages, Clinical Trials, Ethics Review, Participants, Site Management, Visit Schedule, Data Entry (eCRF), Extractions, Adverse Events, Monitoring, Alerts, AI Assistant, Documents, Reports & Exports, Audit Trail, Settings.**
2. Review the dashboard widgets — same layout as other roles, full portfolio scope.

Expected result. All 17 menu entries are present; the dashboard renders with full data.

ACCEPTANCE CHECKLIST

- ☐ The sidebar shows exactly the menu items listed above — no more, no fewer.
- ☐ The welcome header greets you by name.
- ☐ KPI tiles show numbers (not errors or blank tiles).
- ☐ The current date is displayed correctly.
- ☐ All 17 sidebar items are present, including Ethics Review, Monitoring, Audit Trail and Settings.




Fig. 2. Administrator dashboard. Captured from the live system on 12 September 2026.

3 Clinical Trials — full oversight

What this is for. The complete trial portfolio with every lifecycle state. As admin you can drive any transition — including the IEC gate normally reserved for Ethics.

STEPS

1. Click **Clinical Trials** in the sidebar.
2. Confirm all trials render with status chips (draft / iec review / iec approved / ctri registered / active), enrolment bars and the **New Trial** button.
3. Open one trial and confirm the lifecycle card offers the valid next transition.

Expected result. The full portfolio is visible and any trial detail opens with its lifecycle action.

ACCEPTANCE CHECKLIST

- ☐ All trials render with correct status chips.
- ☐ New Trial is available.
- ☐ Trial detail opens with a lifecycle action.

AyuSphere

Research Today · Healthier Tomorrow

Dashboard

Messages

Clinical Trials

Ethics Review

Participants

Site Management

Visit Schedule

Data Entry (eCRF)

Extractions

Adverse Events

Monitoring

Alerts

AI Assistant

Documents

Reports & Exports

Audit Trail

AYUSH

Traditional Wisdom

Modern Evidence

9+

AA

Admin AIAA

Administrator

Search trials, participants, documents...

Clinical Trials

Full lifecycle: protocol → IEC → CTRI → activation → close-out.

New Trial

CTRI-117286

ctri registered

CTRI/2026/09/117286

Comparison of Grape Seed Oil and Aloe Vera Preparation with Clotrimazole for Reducing Oral Candida in Postmenopausal Women with Xerostomia

Grape seed oil with Aloe vera extract · interventional · Phase II · 1 site

ENROLMENT 0/110

0%

CTRI-117131

active

CTRI/2026/09/117131

Comparison of shatapushpa taila matrabasti with tila taila matrabasti in the management of kashtartava or primary dysmenorrhea

Shatapushpa taila; mAtRAbastiH, मात्राबस्ति · interventional · Phase II/III · 1 site

ENROLMENT 14/138

10%

CTRI-116942

iec approved

CTRI/2026/09/116942

Pilot Clinical Evaluation of a Herbo-mineral Hair Removal Cream for Effective Hair Removal

Herbo mineral hair removal cream · interventional · 1 site

ENROLMENT 5/20

25%

CTRI-117332

ctri registered

CTRI/2026/09/117332

A clinical study to evaluate the safety and efficacy of tab Endo-HYLM, a herbomineral formulation in the management of female patients with endometrioma

Endo-HYLM · interventional · Phase I/II · 1 site

ENROLMENT 10/70

14%

CTRI-117085

active

CTRI/2026/09/117085

Role of Panchakarma in treating low back pain

piNDasvedaH, पिण्डस्वेद · interventional · Phase II/III · 1 site

ENROLMENT 0/45

0%

CTRI-116719

ctri registered

CTRI/2026/08/116719

Efficacy of Bala Taila Abhyanga in Occupational Musculoskeletal Disorders among Cobblers: A Pilot Study

Bala Taila; aBya-gga-karma अउंग-कर्म; Yogaraja Guggulu; Maharasnadi Kwatha; Ashwagandha Churna · interventional · 1 site

ENROLMENT 0/15

0%

CTRI-117339

iec approved

CTRI/2026/09/117339

A Study to Compare an Ayurvedic Medicine, Tamra Bhasma with Trikatu Kashaya for Improving Blood Cholesterol Levels

Tamra Bhasam · interventional · Phase IV · 1 site

ENROLMENT 0/80

0%

CTRI-116709

active

CTRI/2026/08/116709

Comparison Of Modified And Classical Takradhara In Management Of Psoriasis Skin Disorder

Vishaghna Mahakashaya; takradhArA, तक्रधरा · interventional · Phase III/IV · 1 site

ENROLMENT 0/24

0%

CTRI-116740

ctri registered

CTRI/2026/08/116740

Role of Panchakarma treatment in Neck pain

nasyam (viShopakramaH), नस्यम् (विशोपक्रम) · interventional · Phase II · 1 site

ENROLMENT 11/40

28%

CTRI-116585

active

CTRI/2026/08/116585

Role of therapeutic Emesis and Purgative treatment on female Hormones and their activity.

vamana-karma, वमन-कर्म; virecana-karma, विरेचन-कर्म · interventional · Phase II · 1 site

ENROLMENT 11/40

28%

Fig. 3. Clinical Trials (admin view) — long page cropped to fit; scroll in the live app to see the rest. Captured from the live system on 12 September 2026.

AyuSphere UAT Guide — confidential demo environment

Page 7 of 14

4 Settings — AI model configuration & user management

What this is for. The admin-only control room. Top card: the AI model that powers Doctor Note extraction and the AI Copilot — provider (Google Gemini or Anthropic Claude), model name and API key; settings saved here override environment variables, and the key is stored server-side and never displayed again. Below: create users with a role, and the account list with activate/deactivate. Every change is audited.

STEPS

1. Click **Settings** in the sidebar — the “Settings — Users & Roles” page opens.
2. In **AI Model**, read the current configuration badge (e.g. *gemini · gemini-2.5-flash · configured in UI*). Note the Provider dropdown (Google Gemini / Anthropic Claude), Model field and API key field (“leave blank to keep current key”). Do NOT save changes unless your test plan provides a spare API key — changing this affects all users' AI features.
3. In **Create user**, fill email, name, role (pi / coordinator / monitor / ethics / pv / admin / regulator) and a password (min 8 chars), then click **Create** to add a test account.
4. In **Accounts**, confirm all 7 demo accounts are listed with role badges, and note the **Deactivate** action per account. Verify you can deactivate (and re-activate) the test account you just created — never the seeded demo accounts.

Expected result. The AI card shows the live provider/model; a new user can be created with any role and immediately log in; deactivation blocks an account's login. All of it lands in the audit trail (settings_ai_update, user create/deactivate).

ACCEPTANCE CHECKLIST

- ☐ AI Model card shows current provider and model with a configured badge.
- ☐ Provider dropdown offers Gemini and Claude; API key is write-only.
- ☐ Create user succeeds with a chosen role.
- ☐ Accounts list shows all users with role badges and Deactivate actions.

AyuSphere

Research Today · Healthier Tomorrow

Dashboard

Messages

Clinical Trials

Ethics Review

Participants

Site Management

Visit Schedule

Data Entry (eCRF)

Extractions

Adverse Events

Monitoring

Alerts

AI Assistant

Documents

Reports & Exports

Audit Trail

AYUSH

Traditional Wisdom
Modern Evidence

Search trials, participants, documents...

9

AA

Admin AIIA
Administrator

Settings — Users & Roles

Account management. Every change is audited.

AI Model

gemini · gemini-2.5-flash · configured in UI

Powers Doctor Note extraction and the AI Copilot. Settings saved here override environment variables; the API key is stored server-side and never displayed again.

PROVIDER

MODEL

API KEY

Google Gemini

gemini-2.5-flash

leave blank to keep current key

Save AI settings

Model examples — Gemini: gemini-2.5-flash, gemini-2.5-pro · Claude: claude-sonnet-5, claude-haiku-4-5-20251001

Create user

EMAIL

NAME

ROLE

PASSWORD (MIN 8)

pi

Create

Accounts (7)

Admin AIIA

admin@aiaa.demo

admin

Deactivate

Ravi Kumar

coordinator@aiaa.demo

coordinator

Deactivate

Prof. S. Iyer

ethics@aiaa.demo

ethics

Deactivate

Meera Nair

monitor@aiaa.demo

monitor

Deactivate

Dr. Ananya Sharma

pi@aiaa.demo

pi

Deactivate

Dr. Farid Khan

pv@aiaa.demo

pv

Deactivate

CDSO Observer

regulator@aiaa.demo

regulator

Deactivate

Fig. 4. Settings — AI Model, create user, and accounts. Captured from the live system on 12 September 2026.

AyuSphere UAT Guide — confidential demo environment

Page 9 of 14

5 Audit Trail

What this is for. The immutable, time-stamped record of every action (ALCOA+): insert-only at the database — nothing here can be edited or deleted. Every row shows when, which actor role, the action (e.g. crf.submit, extraction.approve, trial.transition, message.send, settings.ai_update), the entity, and the before/after change payloads.

STEPS

1. Click **Audit Trail** in the sidebar.
2. Scan recent rows: timestamp, actor role chip, action, entity reference and the BEFORE/AFTER change JSON.
3. Use the filters — **Entity type**, **Entity ID**, **Action** — and click **Filter** to narrow (e.g. action = `extraction.approve`).
4. Perform any small action elsewhere (acknowledge an alert, send a message), return, and confirm it appears at the top of the trail.

Expected result. Every platform action appears with full provenance; filters narrow correctly; there is no edit or delete control anywhere on the page.

ACCEPTANCE CHECKLIST

- ☐ Rows show timestamp, actor role, action, entity and change payload.
- ☐ Filters narrow the list.
- ☐ Your own fresh action appears in the trail.
- ☐ No edit/delete controls exist.

AyuSphere

Research Today · Healthier Tomorrow

Dashboard

Messages

Clinical Trials

Ethics Review

Participants

Site Management

Visit Schedule

Data Entry (eCRF)

Extractions

Adverse Events

Monitoring

Alerts

AI Assistant

Documents

Reports & Exports

Audit Trail

AYUSH

Traditional Wisdom

Modern Evidence

9+

AA

Admin AIIA

Administrator

Q

Search trials, participants, documents...

Audit Trail

Immutable, time-stamped record of every action (ALCOA+). Insert-only at the database — nothing here can be edited or deleted.

ENTITYTYPE

ENTITYID

ACTION

Filter

WHEN	ACTOR ROLE	ACTION	ENTITY	CHANGE
9/12/2026, 3:31:27 PM	admin	message.send	message - ab093c26...	AFTER {"recipientId":"495066d6-d3b3-457d-...
9/12/2026, 3:18:41 PM	pi	message.send	message - 02431728...	AFTER {"recipientId":"92972bdb-0934-4317-...
9/12/2026, 10:56:18 AM	pi	message.send	message - 3feacfe6...	AFTER {"recipientId":"92972bdb-0934-4317-...
9/12/2026, 10:14:07 AM	admin	copilot.query	copilot - copilot...	AFTER {"grounded":true,"question":"Which site...
9/12/2026, 10:13:42 AM	admin	copilot.query	copilot - copilot...	AFTER {"grounded":true,"question":"hello","cit...
9/12/2026, 10:09:17 AM	admin	visit.complete	visit - 24081fae...	BEFORE {"status":"overdue"} AFTER {"status":"completed"}
9/12/2026, 9:59:55 AM	admin	export.dm	trial - d3d235c1...	AFTER {"filename":"CTRI-117253-sdtrm-dm.csv"}
9/12/2026, 9:59:39 AM	admin	export.fhir	trial - fd2fb861...	AFTER {"filename":"CTRI-117133-fhir-bundle.js...
9/12/2026, 9:49:57 AM	admin	extraction.approve	crf_entry - 40c811c8...	AFTER {"source":"extraction","modelId":"gemi...
9/12/2026, 9:49:45 AM	admin	extraction.create	extraction - 83bd39f0...	AFTER {"status":"review","flagCount":1,"fieldCo...
9/12/2026, 9:43:19 AM	admin	extraction.approve	crf_entry - 3dfedc94...	AFTER {"source":"extraction","modelId":"gemi...
9/12/2026, 9:42:32 AM	admin	extraction.create	extraction - 9c6c9ace...	AFTER {"status":"review","flagCount":1,"fieldCo...
9/12/2026, 9:31:56 AM	pi	extraction.approve	crf_entry - 5064dfbf...	AFTER {"source":"extraction","modelId":"gemi...
9/12/2026, 9:30:15 AM	pi	extraction.create	extraction - e354fed3...	AFTER {"status":"review","flagCount":1,"fieldCo...
9/12/2026, 9:23:23 AM	pi	extraction.approve	crf_entry - ea913e12...	AFTER {"source":"extraction","modelId":"gemi...
9/12/2026, 9:23:20 AM	pi	extraction.create	extraction - e504e896...	AFTER {"status":"review","flagCount":0,"fieldC...
9/12/2026, 9:11:21 AM	pi	extraction.create	extraction - fdb3eb39...	AFTER {"status":"review","flagCount":0,"fieldC...
9/12/2026, 8:52:11 AM	pi	extraction.create	extraction - fd99ad02...	AFTER {"status":"review","flagCount":0,"fieldC...
9/12/2026, 8:49:36 AM	pi	extraction.quality_fail	extraction - 67f2c2ed...	AFTER {"score":0.4,"reason":"Image quality","s...
9/12/2026, 8:45:46 AM	pi	extraction.create	extraction - 14ca9c64...	AFTER {"status":"review","flagCount":0,"fieldC...
9/12/2026, 8:41:16 AM	pi	copilot.query	copilot - copilot...	AFTER {"grounded":true,"question":"Summariz...
9/12/2026, 8:40:36 AM	pi	copilot.query	copilot - copilot...	AFTER {"grounded":false,"question":"Any open ...
9/12/2026, 8:14:36 AM	admin	visit.missed	visit - 2da6159c...	BEFORE {"status":"overdue"} AFTER {"status":"missed"}
9/12/2026, 8:14:35 AM	admin	visit.missed	visit - 5bba616a...	BEFORE {"status":"overdue"} AFTER {"status":"missed"}
9/12/2026, 8:14:32 AM	admin	visit.missed	visit - 9e148d4a...	BEFORE {"status":"overdue"} AFTER {"status":"missed"}

Fig. 5. Audit Trail — long page cropped to fit; scroll in the live app to see the rest. Captured from the live system on 12 September 2026.

AyuSphere UAT Guide — confidential demo environment

Page 11 of 14

6 Messages (internal communication)

What this is for. Direct one-to-one messaging with every other user on the platform — coordinators, investigators, monitors, ethics, pharmacovigilance, admin and the regulator — so trial communication stays inside the audited system.

STEPS

1. Click **Messages** in the left sidebar.
2. The **People** list on the left shows every other platform user with their role badge.
3. Click a person (for example the Study Coordinator) to open the conversation panel.
4. Type a short test message (e.g. “UAT test — please ignore”) and send it.
5. Confirm your message appears in the thread with a timestamp.

Expected result. The people list shows the other demo users with role badges. A sent message appears immediately in the conversation and is delivered to the recipient (message events are also recorded in the audit trail).

ACCEPTANCE CHECKLIST

- ☐ Messages page opens from the sidebar.
- ☐ All other demo users are listed with correct role badges.
- ☐ Selecting a person opens a chat panel.
- ☐ A test message can be sent and appears with a timestamp.




Fig. 6. Messages — people list and conversation panel. Captured from the live system on 12 September 2026.

7 Signing out

What this is for. Ends your session so nobody else can act under your account — important on shared computers at trial sites.

STEPS

1. Click your **name / avatar** in the top-right corner of any page.
2. Choose **Sign out** from the menu.

Expected result. You are returned to the login page. Pressing the browser Back button does not re-open protected pages with your data.

ACCEPTANCE CHECKLIST

- ☐ The sign-out option is available from the top-right user menu.
- ☐ After signing out you land on the login page.
- ☐ Navigating directly to a protected page (e.g. /dashboard) redirects you to log in again.

UAT sign-off

Record the outcome of each test case. A test passes when every item of its acceptance checklist is satisfied.

ID	Step / feature	Expected	Pass / Fail	Remarks
TC-01	Logging in	Login succeeds; dashboard opens with correct name and role		
TC-02	Navigation overview & dashboard	Sidebar shows all 17 items; dashboard renders		
TC-03	Clinical Trials — full oversight	Full trial portfolio and lifecycle actions available		
TC-04	Settings — AI model configuration & user management	AI provider/model visible and editable; users can be created and deactivated		
TC-05	Audit Trail	Immutable audit trail lists all actions with before/after; filters work		
TC-06	Messages (internal communication)	People list shows all users; a test message sends and appears in the thread		
TC-07	Signing out	Session ends; protected pages no longer accessible		

Tester name & signature

Date

Reviewed by (QA lead)